

# Intercept® Oral Fluid Drug Testing



## FDA-Cleared, Laboratory-Based Oral Fluid Drug Testing.

Lower Costs with Increased Reliability




### Maximize Productivity and Hiring Time

- Hire employees faster
- Keep employees on the job, minimize lost work time
- Avoid scheduling inconveniences and fees of urine collection
- Eliminate the potential time wasted by shy bladder syndrome and the necessity to maintain a specialized restroom

### Deliver Oral Fluid Convenience

- Speeds up and controls collection process
- Streamline collections done anytime, anywhere
- Preserves employee dignity with non-invasive collection process
- Eliminates adulteration and dilution tactics with directly observed collections every time

### Intercept® Speeds Up the Drug Testing Process



## Intercept® – Laboratory Confidence With Oral Fluid Convenience

	INTERCEPT®	Urine
Days to results for applicant testing <sup>1</sup>	1 day	2-5 days
Time off the job for collection (average) <sup>2</sup>	Under 15 minutes	30-90 minutes
Ecstasy confirmation included in standard panel (NIDA-5) <sup>3</sup>	Yes	No
Prevents tampering, dilution, or masking agents <sup>4</sup>	Yes	No
Certified laboratory analysis	Yes	Yes
Training and certification for collection personnel <sup>5</sup>	<1 hour, CD-based, very easy for existing employees	1-day seminar, many variables on diluted, adulterated or substituted samples

### Swab



**1.** Place pad between cheek and gums for 5 minutes.

### Snap



**2.** Insert pad into open vial, then snap wand.

### Seal



**3.** Recap the vial and apply seal. Send to laboratory with chain of custody form.

For more information call 800.672.7873, or visit us online: [ChooseIntercept.com](http://ChooseIntercept.com)

The Value of a Reliable and Trusted Brand

**Intercept®**

Oral Fluid Drug Test

1. Pre-employment urine drug screen collections and laboratory result reporting times will vary according to individual corporate policies and laboratory arrangements.
2. Oral Fluid Symposium: Advancements and Applications of Oral Fluid Drug Testing; Addiction Professional; Sept, 2004.
3. OraSure Amphetamine-Specific Intercept® MICRO-PLATE EIA Package Insert.
4. Oral fluid drug tests: Effects of adulterants and foodstuffs; Wong, Tran, Tung; Forensic Science International; Volume 150, Issue 2, 10, June 2005; 175-180.
5. Urine collection training times may vary by laboratory vendor.



OraSure Technologies, Inc.

220 East First Street, Bethlehem, PA 18015 USA, phone: 800.ORASURE  
web: [www.OraSure.com](http://www.OraSure.com) • [www.ChooseIntercept.com](http://www.ChooseIntercept.com)

©2008, 2009 OraSure Technologies, Inc. • Intercept® is a registered trademark of OraSure Technologies, Inc.  
Item# INTO109 (Rev. 11/10)

